



## Oxygen Therapy Flow Selector

- **No Calibration required EVER**
- **Lightweight construction**
- **Available with International Fittings**
- **7 Year Warranty**
- **7 Year service Interval**

Excellent performance and robustness, coupled with a very lightweight design.

For use in all Hospital and Emergency environments, the Flow Selector is designed to be plugged into Medical Oxygen Wall Outlets.


The Meditech Flow Selector has been designed to deliver precision flow rates, without the need for any calibration.

Can be supplied in either of two formats: Push-Fit 'Firtree' Fitting for direct attachment on a Oxygen Therapy delivery system With a threaded outlet fitting for the attachment of a humidifier.

### PROBE TYPES:



### THERAPY OUTLET TYPES



**AVAILABLE GAS TYPES:**

- O<sub>2</sub>
- MA4
- O<sub>2</sub>/N<sub>2</sub>O

# SPECIFICATION

Intended Gas: .....	Medical Oxygen, MA4, O2/N2O
Inlet pressure: .....	4 bar (400 kPa)
Set Flow Rates in lpm: .....	Three flow rates options available: 1/2, 1, 2, 3, 4, 6, 8, 10, 15 lpm 1/2, 1, 2, 3, 4, 6, 8, 10, 15, 25 lpm 10, 20, 30, 35, 40, 50, 55, 60, 70 lpm
Inlet Connections: .....	BS Probe + all other international probe types
Outlet (Patient) Connections. ....	Firtree style barbed outlet, 9/16 threaded outlet, 3/8" BSP threaded out- let & 1/4" barb
Intended Life: .....	10 Years
Service Interval: .....	7 Years
Warranty Period: .....	7 Years
Complies With: .....	ISO 15002: 2008

## medITECH MEDICAL GAS EQUIPMENT

**B.N.O.S. Meditech Ltd.**

**9 Fifth Avenue, Bluebridge Industrial Estate,**

**Halstead, Essex CO9 2SZ, England, UK.**

**Tel: +44 (0)1787 479475 • Fax: +44 (0)1787 477747**

**E-mail: [sales@meditech.uk.com](mailto:sales@meditech.uk.com) • [www.meditech.uk.com](http://www.meditech.uk.com)**

**TO PURCHASE OR FOR FURTHER DETAILS:  
Contact us by email or telephone for more details,  
including specific part numbers, prices, etc.**



Meditech Flow Selectors are CE marked in accordance with the European Medical Devices Directive 93/42/EEC and UKCA marked accordance with the UK Medical Devices Regulations 2002. All items in this range are sold with a return to manufacturer warranty against manufacturing defects. These products are part of a series of medical gas and resuscitation products available from B.N.O.S. Meditech. If the product you require is not shown here, or if you require more information, please contact your distributor or contact us directly. Specifications may be subject to change without notice.

B.N.O.S. Meditech Ltd Terms and Conditions of sale apply, available on application or from our website.

Distributor:

**medITECH**  
MARINE

Verlengde Edisonstraat 1,  
8861 NJ Harlingen, Netherlands

E: [sales@meditechmarine.com](mailto:sales@meditechmarine.com)  
T: +31 85 333 22 13

Doc 145M Issue F 29/08/2024