

medITECH

MEDICAL GAS EQUIPMENT

10L Medical Oxygen Cylinder

BNOS Meditech has a Wholesalers Dealers License allowing the sale and distribution of Medical Oxygen. Authorisation Holder Number: WDA (H) 44056

Valve Types:



Pin Index



Bullose



SPECIFICATION

Filling Pressure:.....	200 BAR
Test Pressure:.....	300 BAR
Material:.....	Steel
Diameter:.....	~152mm
Weight(filled with valve):.....	~16Kg
Weight(Empty with valve):.....	~14.7Kg
Height:.....	~ 863mm
Colour Code:.....	White* - EU Standard
Standard:.....	BS EN ISO 7866:2012
Valve:	Pin Index or Bullnose**
Volume (L):.....	2000 Litres
Thread:	25E Tapered

* Other colours available if required

** Other valve types available if required

medITECH

MEDICAL GAS EQUIPMENT

TO PURCHASE OR FOR FURTHER DETAILS:

Contact us by email or telephone for more details,
including specific part numbers, prices, etc.

B.N.O.S. Meditech Ltd.

**9 Fifth Avenue, Bluebridge Industrial Estate,
Halstead, Essex, CO9 2SZ, England, UK.**

TEL: +44 (0) 1787 479 475

FAX: +44 (0) 1787 477747

EMAIL: sales@meditech.uk.com

WEBSITE: www.meditech.uk.com

These products are part of a series of medical gas and resuscitation products available from B.N.O.S. Meditech. If the product you require is not shown here, or if you require more information, please contact us directly. Specifications may be subject to change without notice. B.N.O.S. Meditech Ltd Terms and Conditions of sale apply, available on application or from our website.

Distributor:

medITECH
MARINE

Verlengde Edisonstraat 1,
8861 NJ Harlingen, Netherlands

E: sales@meditechmarine.com

T: +31 85 333 22 13