

# medITECH

MEDICAL GAS EQUIPMENT

## Gas Specific Adaptor

Straight Bar Twin Adaptor



- Available with International Fittings.
- Available with 50mm or 100mm back bar.
- 7 Year Warranty.
- Lightweight construction.
- Service Free.

Excellent performance and robustness, coupled with a very lightweight design.

The Meditech Straight Bar Twin Adaptor is designed for use in all Hospitals and Emergency Enviroments. This device has a Meditech 7 year warranty and is service free for its 10 year recommended life period.

Each Gas Specific Adaptor is available with the probe types and / or the outlet types shown below. If either the probe or outlet type you require is not shown, please contact us to discuss your requirements.

Each Gas Specific Adapter is available with the probe types and / or the outlet types shown below. Other types are available upon request.

### PROBE TYPE



BS PROBE

DISS

DIN

AGA



AFNOR

SIS

CZECH

UNIFOR



### OUTLET TYPES



BS SCHRADER

DISS

DIN

AGA

AFNOR



SIS

CZECH

UNIFOR

OHIO



# SPECIFICATION

|                                   |                                                                                                                                                                                                          |
|-----------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Intended Gas: .....               | <ul style="list-style-type: none"><li>• Medical Oxygen (O2) • Medical Air (MA4)</li><li>• Entonox® (O2/N2O) • Carbon Dioxide (CO2)</li><li>• Surgical Air (MA7) • Nitrous Oxide (N2O) • Vacuum</li></ul> |
| Inlet pressure: .....             | 4 bar (400 kPa)                                                                                                                                                                                          |
| Inlet Connections: .....          | Schrader Probe (BS 5682:2015) and other international connections                                                                                                                                        |
| Outlet Connections (4-BAR): ..... | BS Schrader and Other International Outlets                                                                                                                                                              |
| Intended Life: .....              | 10 Years                                                                                                                                                                                                 |
| Service Interval: .....           | Service free                                                                                                                                                                                             |
| Warranty Period: .....            | 7 Years                                                                                                                                                                                                  |



## TO PURCHASE OR FOR FURTHER DETAILS:

Contact us by email or telephone for more details, including specific part numbers, prices, etc.

**B.N.O.S. Meditech Ltd.**  
**9 Fifth Avenue, Bluebridge Industrial Estate,**  
**Halstead, Essex, CO9 2SZ, England, UK.**  
**TEL: +44 (0) 1787 479 475**  
**FAX: +44 (0) 1787 477747**  
**EMAIL: [sales@meditech.uk.com](mailto:sales@meditech.uk.com)**  
**WEBSITE: [www.meditech.uk.com](http://www.meditech.uk.com)**



Meditech Gas Specific Adaptors are CE marked in accordance with the European Medical Devices Directive 93/42/EEC and UKCA marked accordance with the UK Medical Devices Regulations 2002. All items in this range are sold with a return to manufacturer warranty against manufacturing defects. These products are part of a series of medical gas and resuscitation products available from B.N.O.S. Meditech. If the product you require is not shown here, or if you require more information, please contact your distributor or contact us directly. Specifications may be subject to change without notice.

B.N.O.S. Meditech Ltd Terms and Conditions of sale apply, available on application or from our website.

**medITECH**  
MEDICAL GAS EQUIPMENT

Distributor:

**medITECH**  
MARINE

Verlengde Edisonstraat 1,  
8861 NJ Harlingen, Netherlands

E: [sales@meditechmarine.com](mailto:sales@meditechmarine.com)  
T: +31 85 333 22 13

Doc 796M Issue C 19/08/2024