

# medITECH

MEDICAL GAS EQUIPMENT



## Bobbin Style Flow Selector

- **No Calibration required EVER**
- **Lightweight construction**
- **Available with all International Fittings**
- **2 Year Warranty**
- **10 year service free**
- **Available in Oxygen or Medical Air**
- **Supplied with a fixed probe of your choice or 3/8" BSP coned connector**
- **Supplied with a 9/16" or 3/8" BSP threaded outlet to screw on fir tree connector or humidifier.**

Excellent performance and robustness, coupled with a very lightweight design.

For use in all Hospital and Emergency environments, the Bobbin Style Flow Selector is designed to be plugged into Medical Gas Wall Outlets.

The Meditech Bobbin Style Flow Selector has been designed to deliver precision flow rates, without the need for any calibration.

Can be supplied with 5 flow rates: 0-1lpm, 0-4lpm, 0-15lpm, 0-30lpm & 0-50lpm

Each Bobbin Style Flow Selector is available with the Probe types and Outlet types shown below. If either the Probe type or Outlet type you require is not shown, please contact us to discuss your requirements.

### PROBE TYPE



BS PROBE



DISS



DIN



AGA



AFNOR



SIS



CZECH




UNI



OHMEDA



PURITAN



CHEMETRON



3/8" BSP Nut +  
Liner

### OUTLET TYPES



9/16" THREADED  
THERAPY  
CONNECTOR



3/8" BSP THREADED  
THERAPY  
CONNECTOR

### Available gas Types:

- O2
- MA4

# SPECIFICATION

Intended Gas: .....	Medical Oxygen (O <sub>2</sub> ), Medical Air (MA4)
Inlet pressure: .....	4 bar (400 kPa)
Set Flow Rates in lpm: .....	Three flow rates options available: 0-1lpm, 0-4lpm, 0-15lpm, 0-30lpm & 0-50lpm
Inlet Connections: .....	BS probe or any other international probe types.
Outlet (Patient) Connections. ....	9/16" threaded outlet, 3/8" BSP threaded outlet, supplied with screw on therapy barb
Intended Life: .....	10 Years
Service Interval: .....	10 year service free
Warranty Period: .....	2 Years
Complies With: .....	ISO 15002: 2008



## TO PURCHASE OR FOR FURTHER DETAILS:

Contact us by email or telephone for more details, including specific part numbers, prices, etc.

**B.N.O.S. Meditech Ltd.**  
**9 Fifth Avenue, Bluebridge Industrial Estate,**  
**Halstead, Essex, CO9 2SZ, England, UK.**  
**TEL: +44 (0) 1787 479 475**  
**FAX: +44 (0) 1787 477747**  
**EMAIL: [sales@meditech.uk.com](mailto:sales@meditech.uk.com)**  
**WEBSITE: [www.meditech.uk.com](http://www.meditech.uk.com)**

All items in this range are sold with a return to manufacturer warranty against manufacturing defects. These products are part of a series of medical gas and resuscitation products available from B.N.O.S. Meditech. If the product you require is not shown here, or if you require more information, please contact your distributor or contact us directly. Specifications may be subject to change without notice. B.N.O.S. Meditech Ltd Terms and Conditions of sale apply, available on application or from our website.

# medITECH

MEDICAL GAS EQUIPMENT

Distributor:

**medITECH**  
MARINE

Verlengde Edisonstraat 1,  
8861 NJ Harlingen, Netherlands

E: [sales@meditechmarine.com](mailto:sales@meditechmarine.com)  
T: +31 85 333 22 13

Doc 853M Issue B 14/03/24