

## BEDSIDE MEDICAL OXYGEN SYSTEM

### Oxygen Therapy System (Excluding Cylinder)

**Part Number**  
670-0710-00

With required spec

#### DESCRIPTION

- 1 Oxygen sign
- 2 Twin flow selector on a wall mounting plate (0-15lpm)
- 3 Humidifier x 2
- 4 Adult High Concentration Masks x 10
- 5 40L Medical oxygen cylinder  
(to be purchased seperately)
- 6 Cylinder rack
- 7 O2 Pin Index Regulator, 1 x BS Schrader, Articulating

#### PART NUMBER

- 674-0106-00
- 710-0101-22
- 674-0047-00
- 673-9006-00
- 676-0001-00
- 674-0107-00
- 680-0001-10

**The Oxygen Therapy System (excluding Cylinder)**  
can be tailored to your specifications, contact a  
member of our sales team for more information.

APPROVAL:  
IMO | MFAG



### The Bedside Medical Oxygen System:

In the IMO Bedside solution, the oxygen cylinder is placed between beds, mounted against the bulkhead with the cylinder rack. The regulator with a fixed flow of 50L / min is mounted on the top of the cylinder and the flow selector is mounted on the bulkhead. The double oxygen outlet is fixed to the wall and the hose is connected to the regulator. Two flow selectors, with connections, are placed in the double oxygen wall outlet, which enables the oxygen flow to be regulated between 0-15L/min. A humidifier is connected to each flow selector and this has a disposable oxygen therapy mask attached to it.

# medITECH

## MEDICAL GAS EQUIPMENT

### TO PURCHASE OR FOR FURTHER DETAILS:

**Contact us by email or telephone for more details,  
including specific part numbers, prices, etc.**

B.N.O.S. Meditech Ltd.  
9 Fifth Avenue, Bluebridge Industrial Estate,  
Halstead, Essex CO9 2SZ, England, UK.  
Tel: +44 (0)1787 479475 • Fax: +44 (0)1787 477747  
E-mail: [sales@meditech.uk.com](mailto:sales@meditech.uk.com) • [www.meditech.uk.com](http://www.meditech.uk.com)



Components manufactured by BNOS Meditech and supplied within the kit are CE marked in accordance with the European Medical Devices Directive 93/42/EEC and UKCA marked in accordance with the UK Medical Devices Regulations 2002. All items in this range are sold with a return to manufacturer warranty against manufacturing defects. These products are part of a series of medical gas and resuscitation products available from B.N.O.S. Meditech. If the product you require is not shown here, or if you require more information, please contact your distributor or contact us directly. Specifications may be subject to change without notice.

B.N.O.S. Meditech Ltd Terms and Conditions of sale apply, available on application or from our website.

Distributor:

**medITECH**  
MARINE

E: [sales@meditechmarine.com](mailto:sales@meditechmarine.com)

Doc 854M Issue B 04/09/24