



## microVENT™ INDUSTRIAL RESUSCITATOR

A hand-held, lightweight resuscitator with compact dimensions, the MicroVENT™ Industrial Resuscitator has ergonomic styling to give easy operation with controls in a natural position. It is constructed with a tough reinforced plastic casing to withstand use in difficult situations.

The MicroVENT™ patient valve features 22/15mm connections to take standard airway devices and face masks. The oxygen supply hose is 2 metres long and made of resilient reinforced plastics to retain flexibility in cold conditions. A British Standard (BS) medical oxygen supply probe enables use with all BS standard portable and piped oxygen supplies. Other international hose fittings are also available.

- **Lightweight.**
- **Manual trigger.**
- **Compact.**
- **Suitable for adults only**
- **Easy to use.**
- **Automatic setting preset**
- **Pneumatic Oxygen powered manual Resuscitator.**



# SPECIFICATION

Patient population range: .....	Adult Only
Manual Operation: .....	Automatic & Manual operation. Gas powered.
Automatic setting:.....	Preset @ 600ml - Breath at 12 bpm
Manual flow rate (L/min): .....	40 lpm
Oxygen Concentration V/V: .....	100%
Pressure relief valve with audible warning limits maximum attainable delivery pressure (kPa): .....	4.5 (6.0 on request)
Resuscitator weight excluding supply hose (g):.....	190g
Maximum resuscitator dimensions excluding supply hose (mm):.....	120 x 55 x 100

## TO PURCHASE OR FOR FURTHER DETAILS:

Contact us by email or telephone for more details, including specific part numbers, prices, etc.

### B.N.O.S. Meditech Ltd.

9 Fifth Avenue, Bluebridge Industrial Estate,  
Halstead, Essex, CO9 2SZ, England, UK.

TEL: +44 (0) 1787 479 475

FAX: +44 (0) 1787 477747

EMAIL: [sales@meditech.uk.com](mailto:sales@meditech.uk.com)

WEBSITE: [www.meditech.uk.com](http://www.meditech.uk.com)



Meditech microVENT™s are CE marked in accordance with the European Medical Devices Directive 93/42/EEC and UKCA marked accordance with the UK Medical Devices Regulations 2002. All items in this range are sold with a return to manufacturer warranty against manufacturing defects. These products are part of a series of medical gas and resuscitation products available from B.N.O.S. Meditech. If the product you require is not shown here, or if you require more information, please contact your distributor or contact us directly. Specifications may be subject to change without notice.

B.N.O.S. Meditech Ltd Terms and Conditions of sale apply, available on application or from our website.

**medITECH**  
MEDICAL GAS EQUIPMENT

Distributor:

**medITECH**  
MARINE

Verlengde Edisonstraat 1,  
8861 NJ Harlingen, Netherlands

E: [sales@meditechmarine.com](mailto:sales@meditechmarine.com)  
T: +31 85 333 22 13

Doc 944M Issue C 25/04/2025