

NORWEGIAN SOLUTION MEDICAL OXYGEN SYSTEM

Oxygen Therapy System (Excluding Cylinder)

Part Number
670-0729-00

With required spec

DESCRIPTION

DESCRIPTION	PART NUMBER
1 Oxygen sign	674-0106-00
2 Twin flow selector on a wall mounting plate	710-0101-22
3 Humidifier x 2	674-0047-00
4 Oxygen therapy masks - adult x 10	673-9006-00
5 Cylinder Rack	674-0107-00
6 High pressure Oxygen Hose (1M) Pin Index to UNI 11144/8	131-0463-00
7 40L Medical Oxygen Cylinder (To be purchased separately)	676-0001-00
8 Wall mounting block	131-0464-00
9 Oxygen two stage regulator, UNI 11144/8 Inlet, 1 x BS schrader OUTLET, including pressure relief valve. 0-10 LPM.	674-0204-00

**The Oxygen Therapy System (excluding Cylinder)
can be tailored to your specifications, contact a
member of our sales team for more information.**

Required for vessels with;
• NOR / NIS FLAG
• ISLE OF MAN FLAG



The NS Medical Oxygen System:

In the IMO Norwegian Medical Oxygen Solution, the oxygen cylinder is placed beside the bed and mounted against the bulkhead with the cylinder rack. The regulator is fixed to the wall and connects to the cylinder rack via a high pressure hose assembly. The regulator has a pressure relief valve which can connect onto a pipe which leads to open air allowing any relieved oxygen to be vented out of the room and into the atmosphere.

A twin flow selector is screwed to the bulkhead and plugs into the pressure regulator allowing the user to administer oxygen therapy (via a humidifier if required) to two patients simultaneously at a flow rate between 0 - 15 lpm

10pcs adult oxygen therapy masks are supplied with the kit.

medITECH

MEDICAL GAS EQUIPMENT

TO PURCHASE OR FOR FURTHER DETAILS:
Contact us by email or telephone for more details,
including specific part numbers, prices, etc.

B.N.O.S. Meditech Ltd.
9 Fifth Avenue, Bluebridge Industrial Estate,
Halstead, Essex CO9 2SZ, England, UK.
Tel: +44 (0)1787 479475 • Fax: +44 (0)1787 477747
E-mail: sales@meditech.uk.com • www.meditech.uk.com



Components manufactured by BNOS Meditech and supplied within the kit are CE marked in accordance with the European Medical Devices Directive 93/42/EEC and UKCA marked accordance with the UK Medical Devices Regulations 2002. All items in this range are sold with a return to manufacturer warranty against manufacturing defects. These products are part of a series of medical gas and resuscitation products available from B.N.O.S. Meditech. If the product you require is not shown here, or if you require more information, please contact your distributor or contact us directly. Specifications may be subject to change without notice.

B.N.O.S. Meditech Ltd Terms and Conditions of sale apply, available on application or from our website.

Distributor:

medITECH
MARINE

Verlengde Edisonstraat 1,
8861 NJ Harlingen, Netherlands

E: sales@meditechmarine.com
T: +31 85 333 22 13

Doc 974M Issue B 04/09/24