



Meditech PN:
674-0110-00

Manufacturer PN:
XF-0306G

Autoclavable, Reusable Humidifier

Reusable humidifier for Oxygen Therapy and ventilation / anesthesia inspiratory lines. Non-heated humidifier - ambient temperature functionality with bubble-through humidification system. The units are unbreakable & shatter resistant. Transparent humidification bottle. Graduation shows minimum & maximum water level. Detachable rigid durable polymer cap with gas connectors. Must be capable of disinfection.

SPECIFICATION

Working Volume Level:.....	250ml
Pressure relief safety valve:	>14 kPa (0.1 bar, 2 psi) pressure rating.
Inlet:.....	9/16 th
Outlet:	6mm Barbed Connector
Cap and Connectors:	Made from brass
Bottle and tubes:	Made of polypropylene
Pressure Valve:	Made of brass (chromium plated)
Decontamination procedure:	Autoclave or Clinel wipes
Warranty:	2 years

medITECH

MEDICAL GAS EQUIPMENT

B.N.O.S. Meditech Ltd.

**9 Fifth Avenue, Bluebridge Industrial Estate,
Halstead, Essex CO9 2SZ, England, UK.**

Tel: +44 (0)1787 479475 • Fax: +44 (0)1787 477747

E-mail: sales@meditech.uk.com • www.meditech.uk.com

TO PURCHASE OR FOR FURTHER DETAILS:

Contact us by email or telephone for more details,
including specific part numbers, prices, etc.

Specifications may be subject to change without notice. B.N.O.S. Meditech Ltd. terms
and conditions of sale apply, available on application or from our website.

Distributor:

medITECH
MARINE

E: sales@meditechmarine.com

Doc 979M Issue B 21/10/2025