

Dual Bobbin Style Flow Meter

With BS Probe



- Light weight construction
- Available with International Fittings
- Service free 10 years
- Anodised Aluminium Body

Excellent performance and robustness, coupled with a very lightweight design.

For use in all Hospital and Emergency environments, the Bobbin Style Dual Flow Selector is designed to be plugged into a 4-BAR outlet, such as Medical Oxygen Wall Outlets.

Complies fully with the requirements of ISO15002:2000 Flow Selecting devices for the connection to terminal units of medical gas pipeline systems. Dual Flow Selectors are calibrated to operate from an inlet pressure of 4 BAR (approx. 60psi).

Features;

- Acetal Spindle.
- Polycarbonate inner & outer tube.
- Nitrile O Rings.
- Supplied with plastic tubing nipple.

AVAILABLE GAS TYPES:

- O₂
- MA4

SPECIFICATION

Intended Gas:.....	Medical Oxygen OR Medical Air
Inlet pressure:	4 bar (400 kPa)
Set Flow Rates in lpm:.....	Flow rates options available: 0-1, 0-4, 0-5, 0-15, 0-30 & 0-50lpm
Inlet Connections:	BS Probe (Other international inlets available)
Outlet (Patient) Connections:.....	Firtree style barbed outlet, 9/16 threaded outlet,
Intended Life:.....	10 Years
Service Interval:.....	Service Free
Warranty Period:.....	2 Years
Complies With:.....	ISO 15002: 2000

**TO PURCHASE OR FOR FURTHER DETAILS:
Contact us by email or telephone for more details,
including specific part numbers, prices, etc.**

B.N.O.S. Meditech Ltd.
9 Fifth Avenue, Bluebridge Industrial Estate,
Halstead, Essex CO9 2SZ, England, UK.
Tel: +44 (0)1787 479475 • Fax: +44 (0)1787 477747
E-mail: sales@meditech.uk.com • www.meditech.uk.com



Meditech Flow Selectors are CE marked in accordance with the European Medical Devices Directive 93/42/EEC and UKCA marked in accordance with the UK Medical Devices Regulations 2002. All items in this range are sold with a return to manufacturer warranty against manufacturing defects. These products are part of a series of medical gas and resuscitation products available from B.N.O.S. Meditech. If the product you require is not shown here, or if you require more information, please contact your distributor or contact us directly. Specifications may be subject to change without notice.

B.N.O.S. Meditech Ltd Terms and Conditions of sale apply, available on application or from our website.

medITECH
MEDICAL GAS EQUIPMENT

Distributor:

medITECH
MARINE

Verlengde Edisonstraat 1,
8861 NJ Harlingen, Netherlands

E: sales@meditechmarine.com
T: +31 85 333 22 13

Doc 989M Issue A 17/04/2025